

Veinlite NEO®




Veinlite NEO is a neonatal vein finder for newborns and infants (under 8 weeks), designed for NICU use. It combines through-the-body transillumination with adjustable colour and brightness to help clinicians quickly locate veins and arteries, making difficult IV insertions faster and easier. Its compact design fits comfortably in hand and works well even on the smallest patients.

QUALITY
TECHNOLOGY
AFFORDABLE

Benefits

- Fewer IV access failure
- Less wasted supplies
- Improved patient satisfaction
- Small and compact design



MFG: ISO13485 / ISO9001

NUVEN MEDICA (ASIA PACIFIC) PTE. LTD. SINGAPORE

Office: Zervex Building, No 8, Ubi Road 2, #05-19 Singapore – 408538
www.nuvenmedica.com | michardlau@nuvenmedica.com | +65 9485 0943
TEL: +65 6636 8408/ +65 8940 2900 FAX: +65 6636 8418 EMAIL: sales@nuvenmedica.com

